

OSTEOGEN® PLUG DELIVERY CASE REPORT FOR DOCTORS INFORMATION

Doctors: Please ask your patients if they are allergic to collagen. *Advise patients to avoid alcohol, mouthwash or chlorhexidine for two weeks **as this has been shown to be toxic to fibroblasts** and may retard healing and crestal bridging of soft tissue*.

Fig. 1: Extraction and Debridement

Following anesthesia, extract tooth using standard atraumatic flapless protocol. Thoroughly remove the granulation tissue and the entire pathologic periodontal ligament using the **Impladent Ltd Ultra Coarse Diamond BurSM**, see photo white X, removal of infected ligament. Flush socket and debride twice



See Video on UCD Bur at:

<https://www.impladentltd.com/product-p/bur-kit.htm>



Fig. 2: Generate Bleeding to Establish the RAP

Using the #6 carbide bur included in the UCD kit remove the **Lamina Dura** and make lingual or palatal holes in the lower half of the socket **where trabecular bone is available** to procure medullary blood from the Alveolar Process containing osteoclast cells (220µm) and osteoblasts to trigger the Regional Acceleratory Phenomenon (RAP).⁴ **Profuse bleeding will be absorbed by the hydrophilic OsteoGen® Plug and will help prevent dry socket. Do not hydrate the Plug prior to delivery.**^{5,6}



Fig. 3: Delivery and Initial Plug Compression

Antibiotics with low pH are not conducive to rapid bone formation. **Metronidazole** at pH 6.5 and above (concentration 5mg/mL) is biologically preferred and by injecting into the socket, and mesiodistally at the crest of **interdental septum**.⁷⁻⁸ Hold the OsteoGen® Plug with sterile tweezers, taper Plug apically and deliver into the socket. **Compact the plug aggressively.** The Plug should be large enough, initially with an excess of 3.0mm-5.0mm on average above occlusal plane, so that it can be compressed into and fill the entire socket to the **Soft tissue superior Level**. **Do not place Plug to the crestal bone height!**



Fig. 4: Final Compression: "Making a Membrane"

Plug compression is achieved by using a **Pluggers Instrument to align and compact the bone grafting crystals closer together creating a bioactive membrane barrier which controls migration of connective tissue.**⁹⁻¹¹ Must use more than one Plug for multiple roots. Fill and unite the Plug roots **superiorly at the root trunk to the level of the soft tissue crest.** Leave the top of the Plug intact so it can be compressed into the socket uniformly. No toothbrushing or waterpik for two weeks to yield rapid tissue healing crestally.



Fig. 5: Suturing; Radiolucent to Radiopaque

Passively crisscross suture over Plug, not through the Plug. Do not use Resorbable sutures. Membrane Not Required. Soft tissue should bridge across in 9-14 days*. The non-ceramic OsteoGen® crystals are a low-density graft material. Site will show radiolucent day of placement. Plug resorbs continuously in 3-5 months and longer and is replaced by host bone at a rate depending on the surgery and patient's age or metabolism. The site will become radiopaque and ready for implant placement. OsteoGen® crystals may be used to reinforce the implant osteotomy prior to implant installation.¹¹⁻¹³ After implant installation, place healing screw and place additional crystals crestally to prevent downward migration of epithilium and achieve primary closure.⁹⁻¹¹



Doctor please see videos online at www.impladentltd.com/OsteoGen-Plugs-p/op.htm

BIBLIOGRAPHY

1. Schlesinger C: Extraction with immediate placement and immediate temporization: the ultimate in patient satisfaction. *Implant Practice* 2015; 8(6):20-24
2. Kosinski T: Simple and predictable socket preservation techniques all dentists can implement regardless of extraction or grafting experience. *The Profitable Dentist* Summer 2016:46-50
3. Kosinski T: A Simple and cost-effective socket preservation technique. *Dentistry Today* 2016;35(4)
4. Frost HM: The Regional Acceleratory Phenomenon: a Review. *Henry Ford Hosp Med Journal* 1983; 31(1):3-9
5. Schlesinger C: Using OsteoGen® Strips to graft the gap around an implant: a case study. *Implant Practice* 2017; 10(3):25-30
6. Schlesinger C: A Novel approach grafting around implants. *Dentistry Today* 2017;36(12):50-54
7. Loesche WJ, Grossman N, Giordano J: Metronidazole in periodontitis (IV). The effect of patient compliance on treatment parameters. *J Clin Periodontol.* 1993;20:96-104.
8. Choukroun J, Simonpieri A, Del Corso M, Mazor Z, Sammartino G and Ehrenfest D: Controlling systematic perioperative anaerobic contamination during sinus-lift procedures by using Metronidazole: An innovative approach. *Implant Dentistry* 2008; 17(3):257-270.
9. Spivak JM, Ricci JL, Blumenthal NC, Alexander H. A new canine model to evaluate the biological response of intramedullary bone to implant materials and surfaces. *J Biomed. Mater Research.* 1990; 24:1121-1149.
10. Ricci JL, Blumenthal NC, Spivak JM and Alexander H. Evaluation of a low-temperature calcium phosphate particulate implant material: Physical-chemical properties and in vivo bone response. *J Oral Maxillofacial Surgery.* 1992; 50: 969-978.
11. Valen M and Ganz SD. Part I: A Synthetic Bioactive Resorbable Graft (SBRG) for predictable implant reconstruction. *J Oral Implantology.* 2002; 28(4):167-177.
12. Valen M and Naman LJ: Dental Implant Orthopedics Basic Science Primer with Bone Grafting Techniques for Implant Success (Presented to the AO February 2016).
13. Valen M and Schulman A: Establishment of an implant selection protocol for predetermined success. *J Oral Implantol* 1990;16(3):166-171.